



Recombinant Human Coiled-coil domain-containing protein 60 (CCDC60)

Product Code	CSB-MP818242HU
Abbreviation	CCDC60
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IWA6
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTKVPATKKL QSSPNSGAVR PFYASENLRQ VPDKPMKSIK YMDKEIINLK KDLIRSRFLI QSVKIGRGYF AILREETAKK KKQQQLQKLG EEERNKFQPA EKISEIHYGD TLLSTYDDEK LKTLGARVTR RPFTPIHSCI ISPSLTEAHV EPLFRQLCAL HWLLEALTID HTHHTMKPVI TCWNPKDPGG SKSTIKKINK DKSMGQKWEH FITAPKTKKF KIPTMRVTNR KPSRRGSTLS LSRASGGSSP QSSMISVNPQ SDEPPSVNTQ VTSSKDIEDN ESSSTKPDEE PLYMNLQKLL EMVREDARRT VTIENGMQRK APSILSVLKQ NKSNSAYKEM QTTLKSSERS SSTSAESHIQ PVQKSKNRT NCDINIHYKS GVCNTRAKF YSVAQEAGFC LQDKMEILMK RQEERGIQKF RAFVLVSNFQ KDIAMRHHI SVVKGDAEEI ADHWYFDLLS KLPEDLKNFR PAKKILVKLQ KFGENLDLRI RPHVLLKVLQ DLRIWELCSP DIAVAIEFVR EHHHMPQED YISWLQSRIN IPIGPYSALR
Source	Mammalian cell
Target Names	CCDC60
Protein Names	Recommended name: Coiled-coil domain-containing protein 60
Expression Region	1-550
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.