



Recombinant Human Peptidyl-prolyl cis-trans isomerase-like 6 (PPIL6)

Product Code	CSB-EP818254HU-B
Abbreviation	PPIL6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IXY8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MARPQPCGPP HARGCSPSLP ERPLQVKVVG LFSCPWFQIA KSAEENLKNN HPSKFEDPIL VPLQEFAWHQ YLQEKKRELK NETWEYSSSV ISFVNGQFLG DALDLQKWAH EVWDIVDIKP SALYDALTED FSAKFLRDTK HDFVFLDICI DSSPIGRLIF ELYCDVCPKT CKNFQVLCTG KAGFSQRGIR LHYKNSIFHR IVQNGWIQGG DIVYGKGDNG ESIYGPTFED ENFSVPHNKR GVLGMANKGR HSNGSQFYIT LQATPYLDRK FVAFGQLIEG TEVLKQLELV PTQNERPIHM CRITDSGDPY A
Source	E.coli
Target Names	PPIL6
Protein Names	Recommended name: Peptidyl-prolyl cis-trans isomerase-like 6 Short name= PPIase EC= 5.2.1.8 Alternative name(s): Cyclophilin-like protein PPIL6 Rotamase PPIL6
Expression Region	1-311
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.