



Recombinant Human Procollagen galactosyltransferase 2 (COLGALT2)

Product Code	CSB-YP818260HU
Abbreviation	GLT25D2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8IYK4
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	RFV AERDSEDDGE EPVVFPEspl QSPTVLVAVL ARNAAHTLPH FLGCLERLDY PKSRMAIWAA TDHNVDNTE IFREWLKNVQ RLYHYVEWRP MDEPESYPDE IGPKHWPTSR FAHVMKLRQA ALRTAREKWS DYILFIDVDN FLTNPQTLNL LIAENKTIVA PMLESRLGLYS NFWCGITPKG FYKRTPDYVQ IREWKRTGCF PVPMVHSTFL IDLRKEASDK LTFYPPHQDY TWTFDIIVF AFSSRQAGIQ MYLCNREHYG YLPIPLKPHQ TLQEDIENLI HVQIEAMIDR PPMEPSQYVS VVPKYPDKMG FDEIFMINLK RRKDRRDRML RTLYEQEIEV KIVEAVDGKA LNTSQLKALN IEMLPGYRDP YSSRPLTRGE IGCFLSHYSV WKEVIDRELE KTLVIEDDVR FEHQFKKKLM KLMDNIDQAQ LDWELIYIGR KRMQVKEPEK AVPNVANLVE ADYSYWTLYG VISLEGAQKL VGANPFGKML PVDEFPLVPMY NKHPVAEYKE YYESRDLKAF SAEPLLIYPT HYTGQPGYLS DTETSTIWDN ETVATDWDRT HAWKSRKQSR IYSNAKNTEA LPPPTSLDTV PSRDEL
Source	Yeast
Target Names	COLGALT2
Protein Names	Recommended name: Procollagen galactosyltransferase 2 EC= 2.4.1.50 Alternative name(s): Glycosyltransferase 25 family member 2 Hydroxylysine galactosyltransferase 2
Expression Region	28-626
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.