



# Recombinant Human Centrosomal protein CEP57L1 (CEP57L1)

<b>Product Code</b>	CSB-BP818266HU
<b>Abbreviation</b>	CEP57L1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8IYX8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDSELMHSIV GSYHKPPERV FVPSFTQNEP SQNCHPANLE VTSPKILHSP NSQALILALK TLQEKIHRLE LERTQAEDNL NILSREAAQY KKALENETNE RNLAHQELIK QKKDISIQLS SAQSRCTLLE KQLEYTKRMV LNVEREKNMI LEQQAQLQRE KEQDQMKLYA KLEKLDVLEK ECFRLTTTQK TAEDKIKHLE EKLKEEEHQR KLFQDKASEL QTGLEISKII MSSVSNLKHS KEKKKSSKKT KCIKRRPPWQ ICSKFGALPF VAEKMRQHRD PHILQKPFNV TETRCLPKPS RTTSWCKAIP PDSEKSISIC DNLSSELLMAM QDELQMSME HQELLKQMK TESHSVCDI ECELECLLKK MEIKGEQISK LKKHQDSVCK LQQKVQNSKM SEASGIQQED SYPKGSKNIK NSPRKCLTDT NLFQKNSSFH PIRVHNLQMK LRRDDIMWEQ
<b>Source</b>	Baculovirus
<b>Target Names</b>	CEP57L1
<b>Protein Names</b>	Recommended name: Centrosomal protein CEP57L1 Alternative name(s): Centrosomal protein 57kDa-like protein 1 Centrosomal protein of 57 kDa-related protein Short name= Cep57R Cep57-related protein
<b>Expression Region</b>	1-460
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.