



Recombinant *Saccharomyces bayanus* UPF0357 protein YCL012C (YCL012C)

Product Code	CSB-BP818278SVF
Abbreviation	YCL012C
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8J0M6
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces bayanus</i> (Yeast) (<i>Saccharomyces uvarum</i> x <i>Saccharomyces eubayanus</i>)
Purity	≥85% (SDS-PAGE)
Sequence	NGKVA NKFRELRSRF QSRFNHRHISL NDSFTDDLES GLHSSNFDII SENNNDTRGG LDDISKNEIK QIMEMNNISF DKARLLYMER KFGQNGIAPD GTPIDPKAFT FD
Source	Baculovirus
Target Names	YCL012C
Protein Names	Recommended name: UPF0357 protein YCL012C
Expression Region	26-132
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.