



Recombinant Mouse Amiloride-sensitive amine oxidase [copper-containing] (Abp1)

Product Code	CSB-EP818310MO
Abbreviation	Abp1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8JZQ5
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	<p>VTPPHDKA RIFADLSPQE IKAVHSFLMS RKELGLESSK NLTLAKNSVF LIEMLLPKKK NVLKFLDEGR KSPVREARAI IFFGAQDHPN VTEFAVGPLP RPCYVQALSP RPGHHLSSWSS RPISTAEYDL LYHMLNRAIT PLHQFFLDTT GFSFLGCDDR FLTFTDVAPR GVESGQRRSW LIVQRYVEGY FLHPTGLEIL VDHSSTDVQD WRVEQLWYNG KFYN SPEELA QKYAVGEVEA VVLEEVLED PLPGATEQPP LFSSYKPRGE FHPTVTVAGP HVVQPSGPRY KLEGNVVLYG DWSFSYRLRS SSSLQIFNVL FGGERVAYEV SVQEAVALYG GHTPAGMQTK YIDVGWGLGS VTHELAPGID CPETATFLDA FHYYDSDGPV LYPRALCLFE MPTGVPLRRH FDSNFKGGFN FYAGLKGYVL VLRTTSTVYN YDIWDFIFY PNGVMETKMH ATGYVHATFY TPEGLRHGTR LQTHLLGNIH THLVHYRVDL DVAGTKNSFR TLKTKLENIT NPWSPSHSLV QPTLEQTQYS HEHQAAFRFG QTLPKYLLFS SPQKNRWGHR RSYRLQIHSM AEQVLPPGWQ EERAFTWARY P LAVTKYRES ERYSSSLYNQ NDPWDPPVVF EEFLRNENI ENEDLVAWVT VGFLHIPHSE DVPNTATPGN CVGFLIRPFN FFEEDPSLAS RDTVIVWPQD NGLNHVQRWI PENRDCLVSP PFSYNGTYKP V</p>
Source	E.coli
Target Names	Aoc1
Protein Names	Recommended name: Amiloride-sensitive amine oxidase [copper-containing] Short name= DAO Short name= Diamine oxidase EC= 1.4.3.22 Alternative name(s): Amiloride-binding protein Short name= ABP Histaminase
Expression Region	23-751
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes a membrane glycoprotein that is expressed in many epithelium-rich and/or hematopoietic tissues and oxidatively deaminates



putrescine and histamine. The protein may play a role in controlling the level of histamine and/or putrescine in these tissues. It also binds to and is inhibited by amiloride, a diuretic that acts by closing epithelial sodium ion channels.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.