



# Recombinant Arabidopsis thaliana PP2A regulatory subunit TAP46 (TAP46)

<b>Product Code</b>	CSB-MP818536DOA
<b>Abbreviation</b>	TAP46
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8LDQ4
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGGLAMEEMP LSVLFEQARK IHLAASESGV DQDVVKKGCE MFQKCEDMIG KLALFSSNET KEDISTNNLK YLLVPYYLAE LTKIIQEDR IQIVKASYAK LKEFFSFCEA MELVPDEELE ASSRGGSGAP ADRRALKIAR FKRQKAAEAK LLEIKERKER RGRSTKASAL STPVESGEDD IPDDDSEER EAWLSSINLA ICKAIDLLEM LKREEEMLSA IKERQLKDGE GGFSRDALDD RTKKAETWHR DAAARIQYSK PAQPITCATF AQDVLEGRAS VSQGHEHKNQ PLIFGPASIV GGPLSTERER MIAQVFQPSH RMPTMCIEDA GLTEMNIMND WQEQTKKAIE EATTSWYNDK PLRRKEEDEE DDEDEEAVM KARAFDDWKD DNPRGAGNKK LTPCG
<b>Source</b>	Mammalian cell
<b>Target Names</b>	TAP46
<b>Protein Names</b>	Recommended name: PP2A regulatory subunit TAP46 Alternative name(s): 2A phosphatase-associated protein of 46 kDa
<b>Expression Region</b>	1-405
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.