



Recombinant Arabidopsis thaliana F-box/FBD/LRR-repeat protein At4g00160 (At4g00160)

Product Code	CSB-EP818589DOA-B
Abbreviation	At4g00160
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8LF09
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEQGESLRIR VVMGKDRISE LPDALLIKIL SFLPTKIVVA TSVFSKQWRP LWKLVPNLEF DSEDYDDKEQ YTFSEIVCKS FLSHKAPVLE SFRLEFESEK VDPVDIGLWV GIAFSRHLRE LVLVAADTGT GTAFKFPSSL CTCNTLETLR LVLLILVDIS SPVVMKSLRT LHLELVSYKD ESSIRNLLSG CPILEELLVI RGEDSDIEVF TIDEVPSLKR LTINDDHDGQ EFWGYVINAP SLKYLLIEDL RCPGFCLNAP ELMEANIFDG TSITNDNFLG YLTSVKRLLL NLSPWKITYP TGSIFYQLVS LEMYTREIDW WDLLTLMLEH SPKLQVLKLT DNCVKFHKNG LPGGKWNEPK YVPECLLSHL ETFVWRRFDW GREEEKEIAT YILKNARRLN KATFSTNPVN SEELDKLKER RKVHNELDGV VRASNQCQFV FKFDTSYHVS DSS
Source	E.coli
Target Names	At4g00160
Protein Names	Recommended name: F-box/FBD/LRR-repeat protein At4g00160
Expression Region	1-453
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.