



Recombinant *Oryza sativa* subsp. *japonica* Putative MYST-like histone acetyltransferase 1 (Os07g0626600, LOC_Os07g43360)

Product Code	CSB-EP818597OFG
Abbreviation	Os07g0626600, LOC_Os07g43360
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8LI34
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MGSMEASTAP ENGTAASAAA AASTACNGAG GGGGAAAASN GGGVERRLRS SAASASWASH LPLEVGTRVM CRWRDQKLHP VKVIERRKSS TSSSPADY EY YVHYTEFNRR LDEWVKLEQL DLETVETD VD EKVEDKATSL KMTRHQKRKI DETHVEQGHE ELDAASLREH EEFTKVKNIA KIELGRYEID TWYFSPFPPE YNDSPKLFCC EFCLNFMKRK EQLQRHMKKC DLKHPPGDEI YRSGTLSMFE VDGKKNKVYQ QNLCYLAKLF LDHKTLYYDV DLFLFYVLCE CDDRGCHMVG YFSKEKHSEE SYNLACILTL PPYQRKGYGK FLIAFSYELS KKEGKVGTPPE RPLSDLG LLS YRGYWTRVLL EILKKHKSNI SIKELSDMTA IKADDILSTL QSLDLIQYRK GQHVICADPK VLDRHLKAAG RGGLEVDVSK LIWTPYKEQG
Source	<i>E. coli</i>
Target Names	Os07g0626600
Protein Names	Recommended name: Putative MYST-like histone acetyltransferase 1 EC=2.3.1.48
Expression Region	1-450
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.