



# Recombinant *Oryza sativa* subsp. *japonica* Putative MYST-like histone acetyltransferase 1 (Os07g0626600, LOC\_Os07g43360)

<b>Product Code</b>	CSB-YP818597OFG
<b>Abbreviation</b>	Os07g0626600, LOC_Os07g43360
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8LI34
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MGSMEASTAP ENGTAASAAA AASTACNGAG GGGGAAAASN GGGVERRLRS SAASASWASH LPLEVGTRVM CRWRDQKLHP VKVIERRKSS TSSSPADY EY YVHYTEFNRR LDEWVKLEQL DLETVETD VD EKVEDKATSL KMTRHQKRKI DETHVEQGHE ELDAASLREH EEFTKVKNIA KIELGRYEID TWYFSPFPPE YNDSPKLF FC EFCLNFMKRK EQLQRHMKKC DLKHPPGDEI YRSGTLSMFE VDGKKNKVY G QNLCYLAKLF LDHKTLYYDV DLFLFYVLCE CDDRGCHMVG YFSKEKHSEE SYNLACILTL PPYQRKGYGK FLIAFSYELS KKEGKVGTP E RPLSDLG LLS YRGYWTRVLL EILKKHKSNI SIKELSDMTA IKADDILSTL QSLDLIQYRK GQHVICADPK VLDRHLKAAG RGGLEVDVSK LIWTPYKEQG
<b>Source</b>	Yeast
<b>Target Names</b>	Os07g0626600
<b>Protein Names</b>	Recommended name: Putative MYST-like histone acetyltransferase 1 EC=2.3.1.48
<b>Expression Region</b>	1-450
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.