



Recombinant *Oryza sativa* subsp. *japonica* Putative B3 domain-containing protein Os10g0537100 (Os10g0537100, LOC_Os10g39190)

Product Code	CSB-EP818601OFG
Abbreviation	Os10g0537100, LOC_Os10g39190
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8LNN8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MEFTPISPPT RVAGGEEDSE RGAAAWAVVE KEHMFEEKVVT PSDVGKLNRL VIPKQHAERY FPLDAAAGAG GGGGGGGGGG GKGKGLVLSFE DRTGKAWRFR YSYWNSSQSY VMTKGWSRFV KEKRLGAGDT VSFGRGLGDA ARGRLFIDFR RRRQDAGSFM FPPTAAPP SHHHHQRHHP PLPSVPLCPW RDYTTAYGGG YGYGYGGGST PASSRHVLF LRPQVPAAVL KSVPVHVAAT SAVQEAATTT RPKRVRLFGV NLDCPAAMDD DDDIAGAASR TAASSLLQLP SPSSTSSST AGKKMCSLDL GL
Source	E.coli
Target Names	Os10g0537100
Protein Names	Recommended name: Putative B3 domain-containing protein Os10g0537100
Expression Region	1-312
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.