



# Recombinant Human Tetratricopeptide repeat protein 30B (TTC30B)

<b>Product Code</b>	CSB-EP818703HU
<b>Abbreviation</b>	TTC30B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N4P2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MAGLSGAQIP DGEFTAVVYR LIRNARYAEA VQLLGGELQR SPRSRAGLSL LGYYCYRLQE FALAAECYEQ LGQLHPELEQ YRLYQAQALY KACLYAEATR VAFLLLDNPA YHSRVLRLQA AIKYSEGDLG GSRSLVEQLP SREGGEESGG ENETDGGQINL GCLLYKEGQY EAACSKFFAA LQASGYQPD LSYNLALAYYS SRQYASALKH IAEIIRGIR QHPELGVGMT TEGIDVRSVG NTLVLHQTAL VEAFNLKAAI EYQLRNYEAA QEALDMPPR AEEELDPVTL HNLQALMNMDA RPTGFEKLQ FLLQQNPFPP ETEFGNLLLLY CKYEYFDLAA DVLAENAHLI YKFLTPYLYD FLDAVITCQT APEEAFIKLD GLAGMLTEVL RKLTIQVQEA RHNRDDEAIK KAVNEYDETM EKYIPVMAQ AKIYWNLENY PMVEKIFRKS VEFCNDHDVW KLNVAHVLFM QENKYKEAIG FYEPIVKKHY DNILNVSAIV LANLCVSYIM TSQNEEAEL MRKIEKEEEQ LSYDDPDKKM YHLCIVNLVI GTLYCAKGN YDFGISRVIKS LEPYNKKLGT DTWYYAKRCF LSLLENMSKH TIMLRDSVIQ ECVQFLEHCE LHGRNIPAVI EQPLEEERMH VGKNTVTYES RQLKALIYEI IGWNI
<b>Source</b>	E.coli
<b>Target Names</b>	TTC30B
<b>Protein Names</b>	Recommended name: Tetratricopeptide repeat protein 30B Short name= TPR repeat protein 30B
<b>Expression Region</b>	1-665
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.