



Recombinant Human Tetratricopeptide repeat protein 30B (TTC30B)

Product Code	CSB-EP818703HU
Abbreviation	TTC30B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N4P2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAGLSGAQIP DGEFTAVVYR LIRNARYAEA VQLLGGELQR SPRSRAGLSL LGYYCYRLQE FALAAECYEQ LGQLHPELEQ YRLYQAQALY KACLYAEATR VAFLLLDNPA YHSRVLRLQA AIKYSEGDLP GSRSLVEQLP SREGGEESGG ENETDGQINL GCLLYKEGQY EAACSKFFAA LQASGYQPD LSYNLALAYYS SRQYASALKH IAEIIRGIR QHPELGVGMT TEGIDVRSVG NTLVLHQTAL VEAFNLKAAI EYQLRNYEAA QEALDMPPR AEEELDPVTL HNLALMNMDA RPTGFEKLQ FLLQQNPFPP ETEGNNLLLY CKYEYFDLAA DVLAENAHLI YKFLTPYLYD FLDAVITCQT APEEAFIKLD GLAGMLTEVL RKLTIQVQEA RHNRDDEAIK KAVNEYDETM EKYIPVMAQ AKIYWNLENY PMVEKIFRKS VEFCNDHDVW KLNVAHVLFM QENKYKEAIG FYEPIVKKHY DNILNVAIV LANLCVSYIM TSQNEEAEL MRKIEKEEEEQ LSYDDPDKKM YHLCIVNLVI GTLYCAKGN YDFGISRVIKS LEPYNKLG TDTWYAKRCF LSLENMSKH TIMLRDSVIQ ECVQFLEHCE LHGRNIPAVI EQPLEEERMH VGKNTVTYES RQLKALIYEI IGWNI
Source	E.coli
Target Names	TTC30B
Protein Names	Recommended name: Tetratricopeptide repeat protein 30B Short name= TPR repeat protein 30B
Expression Region	1-665
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.