



# Recombinant Human EF-hand calcium-binding domain-containing protein 4A (EFCAB4A)

<b>Product Code</b>	CSB-EP818705HU-B
<b>Abbreviation</b>	EFCAB4A
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8N4Y2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MASPGKPGAD EAQEEEGELE GGSAGPRAAI LEQAEELFLL CDKEAKGFIT KHDLQGLQSD LPLTPEQLEA VFESLDRAHT GFLTAREFCL GLGMFVGVAS AQGANPC RTP EETFESGGLD VQGTAGSLDE EEEEEERFHT VLEQLGVAPV LGKQRAVRTL WARLQ RERPE LLGSFEDVLI RASACLEEAA RERDGLEQAL RRRESEHERE VRALYEETEQ LREQSRRPPS QNFARGERRS RLELELQSRE QDLERAGLRQ RELEQQLHAQ AAHLEAQAQ NSQLWRAHEA LRTQLEGAQE QIRRLSEAR GRQEQTQRDV VAVSRNMQKE KVSLLRQLEL LRELNTRLRD DRDACEARRA GSSCRKALTT ARLPGPTCCC CCCWARPPRR GSGHLPSAR
<b>Source</b>	E.coli
<b>Target Names</b>	CRACR2B
<b>Protein Names</b>	Recommended name: EF-hand calcium-binding domain-containing protein 4A Alternative name(s): Calcium release-activated calcium channel regulator 2B Short name= CRAC channel regulator 2B
<b>Expression Region</b>	1-399
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.