



Recombinant Human Putative phosphatidylinositol 4-kinase alpha-like protein P1 (PI4KAP1)

Product Code	CSB-MP818731HU
Abbreviation	PI4KAP1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N8J0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MTVEQKFGLF SAEIKEADPL AASEASQPKP CPPEVTPHYI WIDFLVQRFE IAKYCSSDQV EIFSSLLQRS MSLNIGGAKG SMNRHVAAIG PRFKLLTLGL SLLHADVVPN ATIRNVLREK IYSTAFDYFS CPPKFPTQGE KRLREDISIM IKFWTAMFSD KKYLTASQLV PPADIGDLRE QLVEENTGSL SGPAKDFYQR EFDFFNKITN VSAVIKPYPK GDQRKKACLS ALSEVTVQPG CSLPSNPEAI VLDVDYKSGT PM
Source	Mammalian cell
Target Names	PI4KAP1
Protein Names	Recommended name: Putative phosphatidylinositol 4-kinase alpha-like protein P1
Expression Region	1-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.