



Recombinant Human Coiled-coil domain-containing protein 89 (CCDC89)

Product Code	CSB-EP818737HU
Abbreviation	CCDC89
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N998
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MRAPMLQKQQ APRMDTPPPE ERLEKQNEKL NNQEEETEFK ELDGLREALA NLRGLSEER SEKAMLR SRI EEQSQLICIL KRRSDEALER CQILELLNAE LEEKMMQEA E KLKAQGEYSR KLEERFMTLA ANHEMLRFK DEYKSENIKL REENEKLRLE NSSLFSQALK DEEAKVLQLT VRCEALTGEL ETLKERCAQD ACQAQAREKE LLELQSQQAC THTKETEQLR SQLQTLKQQH QQAVEQIAKA EETHSSLSQE LQARLQTVTR EKEELLQLSI ERGKVLQNKQ AEICQLEEKL EIANEDRKHA LERFEQEAVA VDSNLRVREL QRKVDGIQKA YDELRLQSEA FKKHSLDLLS KERELNGKLR HLSP
Source	E.coli
Target Names	CCDC89
Protein Names	Recommended name: Coiled-coil domain-containing protein 89 Alternative name(s): Bc8 orange-interacting protein
Expression Region	1-374
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.