



Recombinant Human Protein FAM71D (FAM71D)

Product Code	CSB-MP818742HU
Abbreviation	FAM71D
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8N9W8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MKKNTSKTTM RINKQDALCT PSHDPRDLQ NMLDGGGEYAP FVSPPMLESN FIQVNRREGES IYLHNRRANWV TVGICFSSST HKIPNVMLLA HLTPGAQKDT ETLFKSLLS PPAEKLVLTR FLPLQFVTLV VHDAENMSLK VKLVSGRAYY LQLCTSAYKQ DTLFSQWVAL ISLLNQEKAK VSKVSEVSSL SGITNSTDIT GSMDVTDVTT FTAILTPYMY AGTGPEHVRD SIDFPEFTDI TDITDVTDLP ENEVPEVPDV RIVTEVIEVR EATEVTDSSD ITNCSGVTVV FENNDLIRAK QEEKEKLKNI LKPGCLQDTK SKSELKESK HVTISNITLT FEGKRYFQTT LTPVESEANT SKEMKDKTSE EKMPDFQSTA LKAEESRSLR TESNTSVLSP HIKSPSNFLK LVPHLSAPFS RE
Source	Mammalian cell
Target Names	FAM71D
Protein Names	Recommended name: Protein FAM71D
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.