



Recombinant Human C2 calcium-dependent domain-containing protein 4A (C2CD4A)

Product Code	CSB-YP818765HU
Abbreviation	C2CD4A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NCU7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MWCLERLRLG PECLRRSGDW LLPGRARGAK SRTTAACANV LTPDRIPEFC IPRLMPRLA LAALRNSWVE EAGMDEGAGR TDWDPRSQAA LSLPHLPRVR TAYGFCALLE SPHTRRKESL LLGGPPAPRP RAHTYGGGGG PDALLGTLRV PRAPGPATPA APGCPRPPQD ALARRPRGCR LLRVPDGLLS RALRAGRSRR LTRVRSVSSG NEDKERRAGS QSPARAPSTS PPSSRVPFPE RLEAEGTVAL GRAGDALRLA AEYCPGTGRL RLRLRAESP AGGAPGPRAV SCRLSLVLRP PGTALRQCST VVGRSRKASF DQDFCFDGLS EDEVRRRLAVR VKARDEGRGR ERGRLLGQGE LSLGALLL
Source	Yeast
Target Names	C2CD4A
Protein Names	Recommended name: C2 calcium-dependent domain-containing protein 4A Alternative name(s): Nuclear-localized factor 1 Protein FAM148A
Expression Region	1-369
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.