



# Recombinant Human Probable aminopeptidase NPEPL1 (NPEPL1)

<b>Product Code</b>	CSB-EP818770HU-B
<b>Abbreviation</b>	NPEPL1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8NDH3
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MANVGLQFQA SAGDSDPQSR PLLLLGQLHH LHRVPWSHVR GKLQPRVTEE LWQAALSTLN PNPTDSCPLY LNYATVAALP CRVSRHNSPS AAHFITRLVR TCLPPGAHRC IVMVCEQPEV FASACALARA FPLFTHRSGA SRRLKKTVT VEFFLVGQDN GPVEVSTLQC LANATDGVRL AARIVDTPCN EMNTDTFLEE INKVGKELGI IPTIIRDEEL KTRGFGGIYG VGKAALHPPA LAVLSHTPDG ATQTIAWVGK GIVYDTGGLS IKGKTTMPGM KRDCGGAAAV LGAFRAAIKQ GFKDNLHAVF CLAENSVGPN ATRPDDIHL L YSGKTVEINN TDAEGRLVLA DGVS YACKDL GADIILDMAT LTGAQGIATG KYHAAVL TNS AEWEAACVKA GRKCGDLVHP LVYCP ELHFS EFTSAVADMK NSVADRDN SP SSCAGLFIAS HIGFDWPGVW VHL DIAAPVH AGERATGFGV ALLLALFGRA SEDPLLNLVS PLGCEVDVEE GDLGRDSKRR RLV
<b>Source</b>	E.coli
<b>Target Names</b>	NPEPL1
<b>Protein Names</b>	Recommended name: Probable aminopeptidase NPEPL1 EC= 3.4.11.- Alternative name(s): Aminopeptidase-like 1
<b>Expression Region</b>	1-523
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.