



Recombinant Human Cytoplasmic polyadenylation element-binding protein 3 (CPEB3)

Product Code	CSB-EP818773HU-B
Abbreviation	CPEB3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NE35
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p> MQDDLMDKS KTQPQPQQQ RQQQPQPES SVSEAPSTPL SSETPKPEEN SAVPALSPAA APPAPNGPDK MQMESPLLPG LSFHQPPQQP PPPQEPAAPG ASLSPSFGST WSTGTTNAVE DSFFQGITPV NGTMLFQNF HHVNPVFGGT FSPQIGLAQT QHHQQPPPPA PAPQPAQPAQ PPQAQPPQQR RSPASPSQAP YAQRSAAAAY GHQPIMTSKP SSSSAVAAAA AAAAASSASS SWNTHQSVNA AWSAPSNPWG GLQAGRDRR AVGVGVGVGV GVPSPNPIS PLKKPFSSNV IAPPKFPRAA PLTSKSWMED NAFRTDNGNN LLPFQDRSRP YDTFNLHSLE NSLMDMIRTD HEPLKKGKHY PSGPPMSFAD IMWRNHFAGR MGINFHHPGT DNIMALNNAF LDDSHGDQAL SSSLSSPTRC QNGERVERYS RKVFGGLPP DIDEDEITAS FRRFGPLVVD WPHKAESKSY FPPKGYAFL FQEESSVQAL IDACLEEDGK LYLCVSSPTI KDKPVQIRPW NLSDSDFVMD GSQPLDPRKT IFVGGVPRPL RAVELAMIMD RLYGGVCYAG IDTDPELKYP KGAGRVAFSN QQSYIAAISA RFVQLQHNDI DKRVEVKPYV LDDQMCDECQ GTRCGGKFAP FFCANVTCLQ YYCEYCWASI HSRAGREFHK PLVKEGGDRP RHVPFRWS </p>
Source	E.coli
Target Names	CPEB3
Protein Names	Recommended name: Cytoplasmic polyadenylation element-binding protein 3 Short name= CPE-BP3 Short name= CPE-binding protein 3 Short name= hCPEB-3
Expression Region	1-698
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.