



Recombinant Human Coiled-coil domain-containing protein 112 (CCDC112)

Product Code	CSB-MP818774HU
Abbreviation	CCDC112
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NEF3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MEKDKHSHFY NQKSDFRIEH SMLEELNKL IHSRKTERAK IQQQLAKIHN NVKKLQHQLK DVKPTPDFVE KLREMMEEIE NAINTFKEEQ RLIYEELIKE EKTTNNELSA ISRKIDTWAL GNSETEKAFR AISSKVPVVK VTPSTLPEEV LDFEKFLQQT GGRQGAWDDY DHQNFVKVRN KHKGKPTFME EVLEHLPGKT QDEVQQHEKW YQKFLALEER KKEIQIWKTK KQKQKREEIF KLKEKADNTP VLFHNKQEDN QKQKEEQRKK QKLAVEAWKK QKSIEMSMKC ASQLKEEEEEK EKKHQKERQR QFKLKLLES YTQQKKEQEE FLRLEKEIRE KAEKAEKRKN AADEISRFQE RDLHKLELKI LDRQAKEDK SQKQRRLAKL KEKVENNVS R DPSRLYKPTK GWEERTKKIG PTGSGPLLHI PHRAIPTWRQ GIQRRV
Source	Mammalian cell
Target Names	CCDC112
Protein Names	Recommended name: Coiled-coil domain-containing protein 112 Alternative name(s): Mutated in bladder cancer protein 1
Expression Region	1-446
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.