



Recombinant Human Protein FAM48A (FAM48A)

Product Code	CSB-EP818775HU
Abbreviation	FAM48A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NEM7
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>MQQALELALD RAEYVIESAR QRPPKRKYLS SGRKSVFQKL YDLYIEECEK EPEVKKLRRN VNLLEKLVQM ETLSCLVVNL YPGNEGYSLM LRGKNGSDSE TIRLPYEEGE LLEYLDAEEL PPILVDLLEK SQVNIFHCGC VIAEIRDYRQ SSNMKSPGYQ SRHILLRPTM QTLICDVHSI TSDNHKWTQE DKLLLESQLI LATAEPLCLD PSIAVTCTAN RLLYNKQKMN TRPMKRCFKR YSRSSLNRQQ DLSHCPPPPQ LRLLDLQKR KERKAGQHYD LKISKAGNCV DMWKRSPCNL AIPSEVDVEK YAKVEKSIKS DDSQPTVWPA HDVKDDYVFE CEAGTQYQKT KLTLQSLGD PLYYGKIQPC KADEESDSQM SPSHSSTDDH SNWFIIGSKT DAERVVNQYQ ELVQNEAKCP VKMSHSSSSGS ASLSQVSPGK ETDQTETVSV QSSVLGKGVK HRPPPIKLPS SSGNSSSSGNY FTPQQTSSFL KSPTPPSSK PSSIPRKSSV DLNQVSMMLSP AALSPASSSQ RTTATQVMAN SAGLNFINVV GSVCGAQALM SGSNPMLGCN TGAITPAGIN LSGLLPSGGL LPNALPSAMQ AASQAGVPFG LKNTSSLRPL NLLQLPGGSL IFNTLQQQQQ QLSQFTPQQP QQPTTCSPQQ PGEQGSEQGS TSQEQALSAQ QAAVINLTGV GSFMQSQAAV LSQLGSAENR PEQSLPQQRF QLSSAFQQQQ QQIQQLRFLQ HQMAMAAAAA QTAQLHHRH TGSQSKSKMK RGTPTTPKF</p>
Source	E.coli
Target Names	SUPT20H
Protein Names	Recommended name: Protein FAM48A Alternative name(s): p38-interacting protein Short name= p38IP
Expression Region	1-779
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.