



Recombinant Human Vacuolar protein sorting-associated protein 37A (VPS37A)

Product Code	CSB-MP818779HU
Abbreviation	VPS37A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NEZ2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MSWLFPLTKS ASSSAAGSPG GLTSLQQQKQ RLIESLRNSH SSIAEIQKDV EYRLPFTINN LTININILLP PQFPQEKPI SVYPPIRHHL MDKQGVYVTS PLVNNFTMHS DLGKIIQSLL DEFWKNPPVL APTSTAFFPYL YSNPSGMSPY ASQGFPLPP YPPQEANRSI TSLSVADTVS SSTTSHTTAK PAAPSGVLS NLPLPIPTVD ASIPTSQNGF GYKMPDVPDA FPELSELSVS QLTDMEQEE VLLEQFLTLP QLKQITDKD DLVKSIEELA RKNLLEPSL EAKRQTVLDK YELLTQMKST FEKKMQRQHE LSESCSASAL QARLKVAAHE AEEESDNIAE DFLEGKMEID DFLSSFMEKR TICHCRRAKE EKLQQAIA MH SQFHAPL
Source	Mammalian cell
Target Names	VPS37A
Protein Names	Recommended name: Vacuolar protein sorting-associated protein 37A Short name= hVps37A Alternative name(s): ESCRT-I complex subunit VPS37A Hepatocellular carcinoma-related protein 1
Expression Region	1-397
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.