



Recombinant Mouse Inositol-trisphosphate 3-kinase A (Itpka)

Product Code	CSB-MP819155MO
Abbreviation	Itpka
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8R071
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MTLPGRPTGM ARPRGAGPCS PGLERAPRRS VGELRLLFEA RCAAIAAAAA AGEPRARGAK RRGQVQPNGL PRAAPAPVIP QLTVTSEEDV TPASPGPPDQ EGNWLPAAGS HLQQPRRLST SSLSTGSSS LLEDSEDDLL SDSESRSRGN VQLETSEDVG QKSHWQKIRT MVNLPVMSPF RKRYSWVQLA GHTGSFKAAG TSGLILKRSS EPEHYCLVRL MADVLRGCVP AFHGIVERDG ESYLQLQDLL DGFDGPCVLD CKMGVRYTYL EELTKARERP KLRKDMYKKM LAVDPEAPTE EEHAQRAVTK PRYMQWREGI SSSTTLGFRI EGIKKADGSC STDFKTTRS R EQVTRVFEEF MQGDAEVLRR YLNRLQQIRD TLEISDFFRR HEVIGSSLLF VHDHCHRAGV WLIDFGKTP LPDGQILDHR RPWEEGNRED GYLLGLDNL I GILASLAER
Source	Mammalian cell
Target Names	Itpka
Protein Names	Recommended name: Inositol-trisphosphate 3-kinase A EC= 2.7.1.127 Alternative name(s): Inositol 1,4,5-trisphosphate 3-kinase A Short name= IP3 3-kinase A Short name= IP3K A Short name= InsP 3-kinase A
Expression Region	1-459
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Target Details	Regulates inositol phosphate metabolism by phosphorylation of second messenger inositol 1,4,5-trisphosphate to Ins(1,3,4,5)P ₄ . The activity of the inositol 1,4,5-trisphosphate 3-kinase is responsible for regulating the levels of a large number of inositol polyphosphates that are important in cellular signaling. Both calcium/calmodulin and protein phosphorylation mechanisms control its activity. It is also a substrate for the cyclic AMP-dependent protein kinase, calcium/calmodulin- dependent protein kinase II, and protein kinase C in vitro.



ITPKA and ITPKB are 68% identical in the C-terminus region.

Reconstitution

We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.