



Recombinant Mouse Ubiquitin-like-conjugating enzyme ATG10 (Atg10)

Product Code	CSB-BP819168MO
Abbreviation	Atg10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8R1P4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MEDEFFGGEKS FQHYCAEFIR HSQQIGDGWE WRTAKECSDG YMCKTQFRIK NEASTPHVGT PASVLTCLPT EENLELPMDD SEVTRPAAVA EVIKHEYHVL YSCSYQVPVL YFRASFLDGR PLALEDIWEG VHECYKPRLL QGPWDTITQQ EHPILGQPFF VLHPCKTNEF MTAVLKNSQK INRNVNYITS WLSLVGPVVG LNLPLSYAKA TSQSE
Source	Baculovirus
Target Names	Atg10
Protein Names	Recommended name: Ubiquitin-like-conjugating enzyme ATG10 EC= 6.3.2.- Alternative name(s): Autophagy-related protein 10 Short name= APG10-like Short name= mAPG10
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.