



Recombinant Mouse SH3 domain-containing kinase-binding protein 1 (Sh3kbp1)

Product Code	CSB-BP819197MO
Abbreviation	Sh3kbp1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8R550
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MVEAIVEFDY QAQHDELTI SVGEVITNIR KEDGGWWEQG INGRRLFPD NFBVREIKKDM KKDLLSNKAP EKPMHDVSSG NALLSSETIL RTNKRGERRR RRCQVAFSYL PQNDDELELK VGDIIIEVVEE VEEGWWEQVL NGKTGMFPSN FIKELSGESD ELGISQDEQL SKSRPEGFLP ASLLPFAHG AKGKTTFFEGT ILYRAAPGKT EGHRRYYSR ETTGSESDGG DSSSTKSEGA NGTMATAAIQ PKKVKGVGFG DIFKDKPIKL RPRSIEVEND FLPVEKTIGK KLPPATSTPD PSKTEMDSRT KTKDYCKVIF PYEAQNDEL TIKEGDIVTL INKDCIDVGW WEGELNGRRG VFPDNFVKLL PSDFDKEGGR PKKPPPSAP VVKQGAGTTE RKHEIKKIPP ERPETLPNRT EEKERPEREP KLDLQKPSVP AIPPKKPRPP KTNSLNRPGA LPPRRPERPV GPLTHTRGDS PKIDLGSAL SGILDKDLSD RSNDIDLEGF DSVISSTEKL SHPTTSRPAK TGRRPPSQSL TSSSLSSPDI FDSPSPEEDK EEHISLAHRG IDVSKKTSKT VTISQVSDNK TSLPPKPGTM AAASSGPASL SSVASSPMSS SLGTAGQRAS SPSLFSTEGK PKMEPAVSSQ AAIEELKMQV RELRTIETM KDQKREIKQ LLSLDEEKK IRLRLQMEVN DIKKALQSK
Source	Baculovirus
Target Names	Sh3kbp1
Protein Names	Recommended name: SH3 domain-containing kinase-binding protein 1 Alternative name(s): Regulator of ubiquitous kinase Short name= Ruk SH3-containing, expressed in tumorigenic astrocytes
Expression Region	1-709
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.