



Recombinant Encephalitozoon cuniculi Cyclin-dependent kinase 1 (CDC28)

Product Code	CSB-MP819379EKH
Abbreviation	CDC28
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8SR86
Product Type	Recombinant Protein
Immunogen Species	Encephalitozoon cuniculi (strain GB-M1) (Microsporidian parasite)
Purity	≥85% (SDS-PAGE)
Sequence	MTESFQKLEK IGEGTYGVVY KARERNTNRV VALKKIRLEN ENEGIPATTI REILLKLNK HSTIVELSDV IYNNNKMYLV FEYVELDLRR YLDRMSDEGR LVEEGFVRKM SQQLLTAMEY CHSRNIFHRD LKPQNILVDP KENIKLADFG LGRAAGIPLR TYTTEVVTLW YRPPELLLCG KYDASVDVW SAACIMAEVV LMRPFPGDS EIDQLFRIFK VLGTPNNSRW SNVENFPNYK VEFVWDPVD LKTIFRGDPD FVDLISKMLE YDPKMRMTAK NGLSHKYFEG MPLIME
Source	Mammalian cell
Target Names	CDC28
Protein Names	Recommended name: Cyclin-dependent kinase 1 Short name= CDK1 Alternative name(s): Cell division control protein 28 homolog EC= 2.7.11.22 Cell division protein kinase 1
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.