



Recombinant Encephalitozoon cuniculi Probable casein kinase I homolog ECU03_0910 (ECU03_0910)

Product Code	CSB-MP819391EKH
Abbreviation	ECU03_0910
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8SS96
Product Type	Recombinant Protein
Immunogen Species	Encephalitozoon cuniculi (strain GB-M1) (Microsporidian parasite)
Purity	>85% (SDS-PAGE)
Sequence	MTTEIRNIKL VQKIASGAFG DIFIGQNTVT NQTVAVKLEK KAHYGQLKHE YGVYKALGGT RTPRIYEYVK ILYENVYVNG LVMELMGKSL EQLFVTCARR FSLKTVLMLG ERMVDNVEYL HHRNYVHRDI KPDNFVFDVQ GDRLYLIDYG LAKEFRNPMT FKHREMRTDK SLTGTARYAS LRTHQGYEQS RRDDLESVGF CMVYFLKGR L PWQGLKAKTK QEKYDRIRES KESISLYELC MGLPKEIHSF CFYVRNLGYE DMPNYAYLRT LLSDALRQRG LRSDGVFDWM VRTPSDSMGD LEIL
Source	Mammalian cell
Target Names	ECU03_0910
Protein Names	Recommended name: Probable casein kinase I homolog ECU03_0910 EC=2.7.11.1
Expression Region	1-304
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.