



Recombinant Human GrpE protein homolog 2, mitochondrial (GRPEL2)

Product Code	CSB-EP819443HU
Abbreviation	GRPEL2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TAA5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	STATQRTA GEDCRSEDPP DELGPPLAER ALRVKAVKLE KEVQDLTVRY QRAIADCENI RRRTQRCVED AKIFGIQSFC KDLVEVADIL EKTTECISEE SEPEDQKLTLEKVFGRLLLLL EAKLKSVFAK HGLEKLTPIG DKYDPHEHEL ICHVPAGVGV QPGTVALVRQ DGYKLHGRTI RLARVEVAVE SQRRRL
Source	E.coli
Target Names	GRPEL2
Protein Names	Recommended name: GrpE protein homolog 2, mitochondrial Alternative name(s): Mt-GrpE#2
Expression Region	33-225
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.