



Recombinant Human Uncharacterized protein C9orf43 (C9orf43)

Product Code	CSB-MP819445HU
Abbreviation	C9orf43
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TAL5
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MDLPDESQWD ETTGGLAVCQ HPQCWATIRR IERGHPRILG SSCKTPLDAE DKLPVLTVVD ILDSGFAAHH LPECTFTKAH SLLSQSSKFY SKFHGRPPKG LPDKSLINCT NRLPKFPVLN LNETQLPCPE DVRNMVVLWI PEETEIHVSQ HGKKKRKNSA VKSKSFLGLS GNQSAGTRVG TPGMIVPPPT PVQLSEQFSS DFLPLWAQSE ALPQDLLKEL LPGGKQTMLC PEMKIKLAMM KKNLPLEKNR PDSVISSKMF LSIHRTLRL PALRYPERLK KLNHLKTEGY RKQQQRQQQQ QQQQKVKTP IKKQEAKKKA KSDPGIQSTS HKHPVTTVHD RLYGYRTLPG QNSDMKQQQQ MEKGTTSKQD STERPKNMYY DHADFHHSVK SPELYETEPT NKDISAPVDA VPEAQAARQK KISFNFSEIM ASTGWNSELK LLRILQDTHDD EDEEDQSSGA E
Source	Mammalian cell
Target Names	C9orf43
Protein Names	Recommended name: Uncharacterized protein C9orf43
Expression Region	1-461
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.