



Recombinant Human Iron-sulfur protein NUBPL (NUBPL)

Product Code	CSB-EP819450HU
Abbreviation	NUBPL
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TB37
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	GA GSETLKQRRT QIMSRGLPKQ KPIEGVKQVI VVASGKGGVG KSTTAVNLAL ALAANDSSKA IGLLDVDVYG PSVPKMMNLK GNPESQSNL MRPLLNYGIA CMSGFLVEE SEPVVWRGLM VMSAIEKLLR QVDWGQLDYL VVDMPPTGD VQLSVSQNIP ITGAVIVSTP QDIALMDAHK GAEMFRRVHV PVLGLVQNMS VFQCPKCKHK THIFGADGAR KLAQTLGLEV LGDIPLHLNI READTGQPI VFSQPESDEA KAYLRIAVEV VRRLPSPSE
Source	E.coli
Target Names	NUBPL
Protein Names	Recommended name: Iron-sulfur protein NUBPL Alternative name(s): IND1 homolog Nucleotide-binding protein-like huInd1
Expression Region	39-319
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.