



Recombinant Human Bromo adjacent homology domain-containing 1 protein (BAHD1)

Product Code	CSB-MP819452HU
Abbreviation	BAHD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8TBE0
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>MTHTRRKSLP MLSSGLTGRR EPLQMEDSNM EQGVEGVEPG MPESPGHLTG RRKNYPLRKR PLVPEKPKAC KVLLTRLENV AGPRSADEAD ELPPDLPKPP SPAPSSDPG LAQPRKRRLA SLNAEALNNL LLEREDTSSL AGTRRSRAGD PHRSRDRDRA TGGWSSSKKR PRLGDLGGGS RDLSPPEPAD EGPRRDGPA PKRLASLNAA AFLKLSQERE LPLRLPRAHA EVDGRSTEPP APKAPRPKWP KVNGKNYPKA WQGASSGEAA GPPGWQGCPD EPWPSATPCG PSVQPSHQPL SKALESPLGL RPHLPLLMGG QAALKPEPGR PGEESPAPKQ ELHQPSFPTP QLSPLPMPGN PADYNGLCVG PELTALGSFY LYCGQEGLQC GGYSPCPMLP EGKLSPVAAAP HEEGLLLAPS SVPSGTPFQH PPWGSSRYCS SEDTGVNGYS ICGVLPLSVT HAGTTCGGCP YKMPFAAEGC RSLGQLEFPL PEAGHPASPA HPLLGCPVPS VPPAAEPVPH LQTPTSEPQT VARACPQSAK PPSGSKSGLR TGSSCRHTAR SKAARRPSHP KQPRVQRPRP RRRRRRRRTNG WVPVGAACEK AVYVLDEPEP AIRKSYQAVE RHGETIRVRD TVLLKSGPRK TSTPYVAKIS ALWENPESGE LMMSLLWYYR PEHLQGGRRSP SMHEPLQNEV FASRHQDQNS VACIEEKCYV LTFAEYCRFC AMAKRREGEL PSRKTALVPP SADYSTPPHR TVPEDTDPEL VFLCRHVYDF RHGRILKNPQ</p>
Source	Mammalian cell
Target Names	BAHD1
Protein Names	Recommended name: Bromo adjacent homology domain-containing 1 protein Short name= BAH domain-containing protein 1
Expression Region	1-780
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.