



# Recombinant Human Phosphatidylinositol 5-phosphate 4-kinase type-2 gamma (PIP4K2C)

<b>Product Code</b>	CSB-MP819455HU
<b>Abbreviation</b>	PIP4K2C
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8TBX8
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	ASSSVPPAT VSAATAGPGP GFGFASKTKK KHFVQQKVKV FRAADPLVGV FLWGVVVAHSIN ELSQVPPPVM LLPDDFKASS KIKVNNHLFH RENLPSHFKE KEYCPQVFRN LRDRFGIDDQ DYLVSLTRNP PSESEGSDGR FLISYDRTL I KEVSSEDIA DMHSNLSNYH QYIVKCHGNT LLPQFLGMYR VSVDNEDSYM LVMRNMFSHR LPVHRKYDLK GSLVSREASD KEVKELPTL KDMDFLNKNQ KVYIGEEEEK IFLEKLKRDV EFLVQLKIMD YSLLLGIHDI IRGSEPEEEA PVREDESEVD GDCSLTGPPA LVGSYGT SPE GIGGYIHSR PLGPGFESF IDVYAIRSAE GAPQKEVYFM GLIDILTQYD AKKKAHAHA TVKHGAGAEI STVHPEQYAK RFLDFITNIF A
<b>Source</b>	Mammalian cell
<b>Target Names</b>	PIP4K2C
<b>Protein Names</b>	Recommended name: Phosphatidylinositol 5-phosphate 4-kinase type-2 gamma EC= 2.7.1.149 Alternative name(s): Phosphatidylinositol 5-phosphate 4-kinase type II gamma Short name= PI(5)P 4-kinase type II gamma Short name= PIP4KII-
<b>Expression Region</b>	2-421
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.