



# Recombinant Human Mirror-image polydactyly gene 1 protein (MIPOL1)

<b>Product Code</b>	CSB-MP819463HU
<b>Abbreviation</b>	MIPOL1
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8TD10
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MENWSKDITH SYLEQETTGI NKSTQPDEQL TMNSEKSMHR KSTELVNEIT CENTEWPGQR STNFQIISSY PDDESUYCTT EKYNVMEHRH NDMHYECMTP CQVTSDDSDKE KTIAFLLKEL DILRTSNKKL QQKLAKEDKE QRKLKFKLEL QEKETEAKIA EKTAALVEEV YFAQKERDEA VMSRLQLAIE ERDEAIARAK HMESLKVLE NINPEENDMT LQELLNRINN ADTGIAIQKN GAIIVDRIYK TKECKMRITA EEMSALIEER DAALSKCKRL EQELHHVKEQ NQTSANNMRH LTAENQERA LKAKLLSMQQ ARETAVQQYK KLEEEIQT LR VYYSLHKLSL QEENLKDQFN YTLSTYEEAL KNRENIVSIT QQQNEELATQ LQQALTERAN MELQLQHARE ASQVANKEVQ KLERLVDVLR KKVGTGMT RT VI
<b>Source</b>	Mammalian cell
<b>Target Names</b>	MIPOL1
<b>Protein Names</b>	Recommended name: Mirror-image polydactyly gene 1 protein
<b>Expression Region</b>	1-442
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.