



Recombinant Mouse Notchless protein homolog 1 (Nle1)

Product Code	CSB-EP819766MO-B
Abbreviation	Nle1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VEJ4
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	AAAVVEEAA AGDVQRLLVQ FQDEGGQLLG SPFDVPVDIT PDKLQLVCNA LLAQEEPLPL AFYVHDAEIV SSLGKTLESQ SVETEKIVDI IYQPQAVFRV RAVTRCTSSL EGHSEAVISV AFSPTGKYLA SGSGDTTVRF WDLSTETPHF TCKGHRHWVL SISWSPDGKK LASGCKNGQV LLWDPSTGLQ VGRTLTHGHSK WITGLSWEPL HMNPECRYVA SSSKDGSVRV WDTTAGRCER ILTGHTQSVT CLRWGGDGLL YSASQDRTIK VWRAHDGVLC RTLQGHGHWW NTMALSTDYA LRTGAFEPAE ATVNAQDLQG SLKELKERAS SRYNLVRGQG PERLVSGSDD FTLFLWSPA E DKKPLARMTG HQALINQVLF SPDSRIVASA SFDKSIKLWD GRTGKYLASL RGHVAAVYQI AWSADSRLLV SGSSDSTLKV WDVKAQKLAT DLPGHAEVY AVDWSPDGQR VASGGKDKCL RIWRR
Source	E.coli
Target Names	Nle1
Protein Names	Recommended name: Notchless protein homolog 1
Expression Region	2-485
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.