



Recombinant Arabidopsis thaliana Zinc finger AN1 and C2H2 domain-containing stress-associated protein 11 (SAP11)

Product Code	CSB-BP819816DOA
Abbreviation	SAP11
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8VZ42
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MGTPEFPDLG KHCSVDVCKQ IDFLPFTCDR CLQVFCLDHR SYMKHSCP KG DREDVTVIC PLCAKGVRLN PNEDPNITWE KHVNTDCDPS NYEKATK KKK CPVPRCKEYL TFSNTIKCRD CNVDHCLKHR FGDHTCPGP RKLPMGF LS SSTTRKEAKT TRPNKAHPST SSSSSSRWS NLLSSAEAGI SRLGNDISQK LQFSSSKDNG IVEVCPQCGA KFSSVTSLVE HVEK THERNK KQNHGNVTVD VCPRCSRGRF DPVDLVNHIE RDHRGTSKA
Source	Baculovirus
Target Names	SAP11
Protein Names	Recommended name: Zinc finger AN1 and C2H2 domain-containing stress-associated protein 11 Short name= AtSAP11
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.