



Recombinant Arabidopsis thaliana

Uncharacterized protein At4g06744 (At4g06744)

Product Code	CSB-YP819831DOA
Abbreviation	At4g06744
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8W3M4
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	LPDR GTQQLTNDGL TGARREILSF LDQRLAVVYP VIQRFKSLIT LDPYNVTKTW IGSDICSYRG FHCDNPPHNK TAVTVASIDF NGFQLSAPSI EGFIDQFADL ALFHVNSNNF GGTVPSKIVN LRYLYELDIS NNRFTGQFPT AVVGMSGLTF IDIRFNSFSG SIPPQILGQN LEVLFINNDG FTASLPEIPG DGTTHILFLT LANNKFNGPL PRSILRSMST LTEVLFLNND FTGCIPHEIG FLTGASVIDI GGNKLTGPLP LSLMCLEKVE QLNFAFNLLF GAVPEAVCML LRDNLVNLSL SDNYFTHVGP WCRGELLEKGV LDVGNNCIPF FPGQRSMQEC AEFFVKPKKY YCPHMWFHSF FSCRYSHISS SSSSAFMPMV APSP
Source	Yeast
Target Names	At4g06744
Protein Names	Recommended name: Uncharacterized protein At4g06744
Expression Region	17-404
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.