



# Recombinant Arabidopsis thaliana Probable calcium-binding protein CML49 (CML49)

<b>Product Code</b>	CSB-MP819837DOA
<b>Abbreviation</b>	CML49
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8W4L0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Arabidopsis thaliana (Mouse-ear cress)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSGYPPSSQG YGYGGNPPPP QPPYGSTGNN PPPYGSSGSN PPPPYGSSAS SPYAVPYGAQ PAPYGAPPSA PYASLPGDHN KPHKEKPHGA SYGSPSPGGY GAHPSSGSPD YGGYGGAPQQ SGHGGGYGGA PQQSGHGGGY GAPPQASYG SPFASLVPSA FPPGTDPNIV ACFQAADRDN SGFIDDKELQ GALSSYNQSF SIRTVHLLMY LFTNSNVRKI GPKEFTSLFF SLQNWRSIFE RFDKDRSGRI DTNELRDALM SLGFSVSPVI LDLLVSKFDK SGGRNRAIEY DNFIECCLTV KGLTEKFKEK DTALSGSAIF NYENFMLTVL PFLVA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	CML49
<b>Protein Names</b>	Recommended name: Probable calcium-binding protein CML49 Alternative name(s): Calmodulin-like protein 49
<b>Expression Region</b>	1-335
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.