



Recombinant GDP-fucose protein O-fucosyltransferase 2 (pad-2)

Product Code	CSB-YP819879CXY
Abbreviation	pad-2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8WR51
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	ENSDQTVSRV DSNRYSVAAE KKFLLYDVNF GEGFNLRRDV YMRVANTVRS LRDSGENYIL VLPPWGRLLH WKRMEVALSW RLFFDLESLN RFIPVIEFED FLDENRPIDQ VIYLQHYAEG WGTEYVRKFE KRSCLPPAES HYKQVEEFKW KGFYFYSYEDV YSRNFQCVSI QGDSGTLKDL LKHSNFSEST SIMVDRAETI LHEHYGEVDY WKARRSMRYS NDLVDVADAF RKKYLDSDDK RDKTKLVDDW TKEKPRRTAI GGPYLGIIHWR RRDFLYARRA QLPTIPGTAK ILQDLCKKLD LQKIYLATDA PDQEVDELKA LLNGELEVYR FTDTQKLNDG QIAIIDQYLC AHAAYFIGSY ESTFTFRIQE DREIIGFPIS TTFNRLCPDT EPTCEQPAKW KIVY
Source	Yeast
Target Names	pad-2
Protein Names	Recommended name: GDP-fucose protein O-fucosyltransferase 2 EC= 2.4.1.221 Alternative name(s): Patterning defective protein 2 Peptide-O-fucosyltransferase 2 Short name= O-FucT-2
Expression Region	21-424
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.