



# Recombinant Multidrug resistance outer membrane protein MdtQ (mdtQ)

<b>Product Code</b>	CSB-EP820214EOD
<b>Abbreviation</b>	mdtQ
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q8X659
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O157:H7
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	CAPMHETRQ ALSQQTPAAQ VDTALPTALK NGWPDSQWWL EYHDNQLTSL INNALQNAPD MQVAEQRIQL AEAQAKAVAT QDGPQIDFSA DMERQKMSAE GLMGPFALND PAAGTTGPWY TNGTFGLTAG WHLDIWGKNR AEVTARLGTV KARAAEREQT RQLLAGSVAR LYWEWQTQAA LNTVLQQIEK EQNTIIATDR QLYQNGITSS VEGVETDINA SKTRQQLNDV AGKMKIIEAR LIALTNHQTK SLTLKPVALP KVASQLPDEL GYSLLARRAD LQAAHWYVES SLSTIDAAKA AFYDPINLMA FLQQDALHLS DLFRRHSAQQM GVTAGLTLPI FDSGRLNANL DIAKAESNLS IASYNKAVVE AVNDVARAAS QVQTLAEKNQ HQAQIERDAL RVVGLAQARF NAGIAGSRV SEARIPALRE RANGLLLQGQ WLDASIQLTG ALGGGYKR
<b>Source</b>	E.coli
<b>Target Names</b>	mdtQ
<b>Protein Names</b>	Recommended name: Multidrug resistance outer membrane protein MdtQ
<b>Expression Region</b>	22-478
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.