



Recombinant HTH-type transcriptional regulator hdfR (hdfR)

Product Code	CSB-BP820257EOD
Abbreviation	hdfR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8XAW1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O157:H7
Purity	≥85% (SDS-PAGE)
Sequence	MDTELLKTFL EVSRTRHFGR AAESLYLTQS AVSFRI RQLE NQLGVNLFTR HRNNIRLTAA GEKLLPYAET LMSTWQAARK EVAHTSRHNE FSIGASASLW ECMLNQWLGR LYQNQDAHTG LQFEARIAQR QSLVKQLHER QLDLLITTEA PKMDEFSSQL LGYFTLALYT SAPSKLKGDL NYLRLEWGPD FQQHEAGLIG ADEVPILTTS SAELAQQQIA MLNGCTWLPV SWARKKGGLH TVVDSTTLRS PLYAIWLQNS DKNTLIRDLL KINVLDEVY
Source	Baculovirus
Target Names	hdfR
Protein Names	Recommended name: HTH-type transcriptional regulator hdfR Alternative name(s): H-NS-dependent flhDC regulator
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.