



Recombinant Mouse Inorganic pyrophosphatase 2, mitochondrial (Ppa2)

Product Code	CSB-EP820989MO
Abbreviation	Ppa2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q91VM9
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	SLY RTEELGHPRS QDYRLFFKHV AGHYISPFHD IPLKADCKEE HDIPRKKARN DEYENLFNMV VEIPRWTNAK MEIATEEPLN PIKQDIKNGK LRYTPNIFPH KGYIWNYGAL PQTWEDPHLR DKSTDCCGDN DPIDVCEIGS KVL SRGDVVH VKILGTLALI DQSETDWKII AINVNDPEAE KFHDIDDVKK FKPGYLEATL NWFRLYKVPD GKPENKFAFN GEFKNKAFAL DVINSAHERW KEMVMKKCDK GAISCVNVHI CDSPFHCTME EARSLVESVP TPSMNKESNV EEEVWHFLRN
Source	E.coli
Target Names	Ppa2
Protein Names	Recommended name: Inorganic pyrophosphatase 2, mitochondrial EC= 3.6.1.1 Alternative name(s): Pyrophosphate phospho-hydrolase 2 Short name= PPase 2
Expression Region	28-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This protein is localized to the mitochondrion, is highly similar to members of the inorganic pyrophosphatase (PPase) family, and contains the signature sequence essential for the catalytic activity of PPase. PPases catalyze the hydrolysis of pyrophosphate to inorganic phosphate, which is important for the phosphate metabolism of cells. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.