



# Recombinant Schizosaccharomyces pombe Ubiquitin carboxyl-terminal hydrolase 6 (ubp6)

<b>Product Code</b>	CSB-MP821650SXV
<b>Abbreviation</b>	ubp6
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q92353
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MMIPIAIRWQ GKKYDLEIEP NETGSTLKHQ LYSLTQVPPE RQKVIVKGGQ LKDDVLLGSV GIKPNATLLM MGTAGELPTA MPIPAVESVE QEESEDDGYP SGLINLGNTC YMNSTVQMLR APELSDAVS QFNSSGGLVA EYRTLLNSMQ SNAPVTPMRF LQSLRMEYPQ FAEMSRETGG YAQQDAEECW SFLLSVLQRS LSSEWVQKNM AGKLLSTMKC DENEVQEQPS ISHDTFLSLP CHISMHTSYM TQGILEGLTQ KISKHSDVLN RDAMYSKISR ISRLPNYLTV NFVRFYWKAS IGKKAKILRK VKFPFELDAV EFCTPELSQK LIPVRDKLRE IEKNDEEHER AAKRIKIQPS EDEKEAEAEC RLTQVATCQS LVDPELADDE GANPTGLYDL VGVLSHAGAS ASSGHYQAWI RNSNNRAEFW RFNDAKVSIV PAEKIETLDG GGEADSAYIL LYKAKDIA
<b>Source</b>	Mammalian cell
<b>Target Names</b>	ubp6
<b>Protein Names</b>	Recommended name: Ubiquitin carboxyl-terminal hydrolase 6 EC= 3.4.19.12 Alternative name(s): Deubiquitinating enzyme 6 Ubiquitin thiolesterase 6 Ubiquitin-specific-processing protease 6
<b>Expression Region</b>	1-468
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.