



Recombinant Human RNA polymerase-associated protein RTF1 homolog (RTF1)

Product Code	CSB-EP821675HU-B
Abbreviation	RTF1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92541
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MRGRLCVGRA AAAAAAVAVP LAGGQEGSPG GGRRGSRGTT MVKKRKGRVV IDSDTEDSGS DENLDQELLS LAKRKRSDSE EKEPPVSQPA ASSDSETS DS DDEWTFGSNK NKKKGKARKI EKKGTMKKQA NKTASSGSSD KDSSAESSAP EEGEVSDSDS NSSSSSSSDSD SSEDEEFHD GYGEDLMGDE EDRARLEQMT EKEREQELFN RIEKREVLKR RFEIKKLLKT AKKKEKKEKK KKQEEQEKK KLTQIQESQV TSHNKERRSK RDEKLDKKSQ AMEELKAERE KRKNRTAELL AKKQPLKTSE VYSDDEEEEE DDKSSEKSDR SSRTSSSDEE EEKEEIPPKS QPVS LPEELN RVRLSRHLE RWCHMPFFAK TVTGCFVRIG IGNHNSKPVY RVAEITGVVE TAKVYQLGGT RTNKGLQLRH GNDQRVFRLE FVSNQEFTES EFMKWKEAMF SAGMQLPTLD EINKKELSIK EALNYKFNDQ DIEEIVKEKE RFRKAPPNYA MKKTQLLKEK AMAEDLGDQD KAKQIQDQLN ELEERAEALD RQRKNISAI SYINQRNREW NIVESEKALV AESHNMKNQQ MDPFTRRQCK PTIVSNSRDP AVQAAILAQL NAKYGSGVLP DAPKEMSKGQ GKDKDLNSKS ASDLSEDLFK VHDFDVKIDL QVPSSSEKAL AITSKAPPAK DGAPRRSLNL EDYKKRRGLI
Source	E.coli
Target Names	RTF1
Protein Names	Recommended name: RNA polymerase-associated protein RTF1 homolog
Expression Region	1-710
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final



concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.