



# Recombinant Mouse Interleukin-24 (IL24)

|                          |  |
|--------------------------|--|
| <b>Product Code</b>      | CSB-BP821683MO   |
| <b>Abbreviation</b>      | IL24   |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q925S4   |
| <b>Product Type</b>      | Recombinant Protein  |
| <b>Immunogen Species</b> | Mus musculus (Mouse)   |
| <b>Purity</b>            | ≥85% (SDS-PAGE)  |
| <b>Sequence</b>          | QEFR SGSCQVTGVV LPELWEAFWT VKNTVQTQDD ITSIRLLKPQ<br>VLRNVSGAES CYLAHSLKLF YLNTVFKNYH SKIAKFKVLR SFSTLANNFI<br>VIMSQLQPSK DNSMLPISES AHQRFLFRR AFKQLDTEVA LVKAFGEVDI<br>LLTWMQKFYH L  |
| <b>Source</b>            | Baculovirus  |
| <b>Target Names</b>      | IL24   |
| <b>Protein Names</b>     | Recommended name: Interleukin-24 Short name= IL-24 Alternative name(s):<br>IL-4-induced secreted protein Melanoma differentiation-associated gene 7<br>protein Short name= MDA-7 Th2-specific cytokine FISP  |
| <b>Expression Region</b> | 27-181   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.  |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.  |
| <b>Protein Length</b>    | Full Length of Mature Protein  |
| <b>Target Details</b>    | This gene encodes a member of the IL10 family of cytokines. It was identified as a gene induced during terminal differentiation in melanoma cells. This protein can induce apoptosis selectively in various cancer cells. Overexpression of this gene leads to elevated expression of several GADD family genes, which correlates with the induction of apoptosis. The phosphorylation of mitogen-activated protein kinase 14 (MAPK7/P38), and heat shock 27kDa protein 1 (HSPB2/HSP27) are found to be induced by this gene in melanoma cells, but not in normal immortal melanocytes. Alternatively spliced transcript variants encoding distinct isoforms have been reported. |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.  |



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.