



# Recombinant Human La-related protein 4B (LARP4B)

<b>Product Code</b>	CSB-EP821686HU-B
<b>Abbreviation</b>	LARP4B
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q92615
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MTSDQDAKVV AEPQTQRVQE GKDSAHLMNG PISQTTSQTS SIPPLSQVPA TKVSELNPNA EVWGAPVLHL EASSAADGVS AAWEEVAGHH ADRGPQGSDA NGDGDQGHEN AALPDPQESD PADMNALALG PSEYDSL PEN SETGGNESQP DSQEDPREVL KKTLEFCLSR ENLASDMYLI SQMDSQYVP ITTVANLDHI KKLSTDVDLI VEVLRSLPLV QVDEKGEKVR PNQNR CIVIL REISESTPVE EVEALFKGDN LPKFINCEFA YNDNWFITFE TEADAQQAYK YLREEVKTFQ GKPIKARIKA KAI AINTFLP KNGFRPLDVS LYAQQRYATS FYFPPMYS PQ QFPLYSLIT PQTWSATHSY LDPPLVTPFP NTGFINGFTS PAFKPAASPL TSLRQYPPRS RNPSKSHLRH AIPSAERGPG LLES PSIFNF TADRLINGVR SPQTRQAGQT RTRIQNPSAY AKREAGPGRV EPGSLESSPG LGRGRKNSFG YRKKREEKFT SSQTQSPTPP KPPSPSFELG LSSFPLPGA AGNLKTEDLF ENRLS SLIIG PSKERTLSAD ASVNTLPVVV SREPSVPASC AVSATYERSP SPAHL PDDPK VAEKQRETHS VDRLPSALTA TACKSVQVNG AATELRKPSY AEICQRTSKE PPSSPLQPQK EQKPNTVGCG KEEKLAEPA ERYREPPALK STPGAPRDQR RPAGGRPSPS AMGKRLSREQ STPPKSPQ
<b>Source</b>	E.coli
<b>Target Names</b>	LARP4B
<b>Protein Names</b>	Recommended name: La-related protein 4B Alternative name(s): La ribonucleoprotein domain family member 4B La ribonucleoprotein domain family member 5 La-related protein 5
<b>Expression Region</b>	1-738
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the



contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

### Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.