



Recombinant Human snRNA-activating protein complex subunit 3 (SNAPC3)

Product Code	CSB-MP821717HU
Abbreviation	SNAPC3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92966
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAEGSRGGPT CSGVGGRQDP VSGSGGCNFP EYELPELNTR AFHVGAFGEL WRGRLRGAGD LSLREPPASA LPGSQAADSD REDAAVARDL DCSLEAAAE L RAVCGLDKLK CLEDGEDPEV IPENTDLVTL GVRKRFLEHR EETITIDRAC RQETFVYEME SHAIGKKPEN SADMIEEGEL ILSVNILYPV IFHKHKEHKP YQTMLVLGSQ KLTQLRDSIR CVSDLQIGGE FSNTPDQAPE HISKDLYKSA FFYFEGTFYN DKRYPECRDL SRTIIEWSSES HDRGYGKFQT ARMEDFTFND LCIKLGFPYL YCHQGDCEHV IVITDIRLVH HDDCLDRTLY PLLIKKHWLW TRKCFVCKMY TARWVTNND S FAPEDPCFFC DVCFRMLHYD SEGNKLGEFL AYPYVDPGTF N
Source	Mammalian cell
Target Names	SNAPC3
Protein Names	Recommended name: snRNA-activating protein complex subunit 3 Short name= SNAPc subunit 3 Alternative name(s): Proximal sequence element-binding transcription factor subunit beta Short name= PSE-binding factor subunit beta Short n
Expression Region	1-411
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.