



Recombinant Rickettsia conorii Succinate dehydrogenase iron-sulfur subunit (sdhB)

Product Code	CSB-MP821823RMS
Abbreviation	sdhB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q92JJ8
Product Type	Recombinant Protein
Immunogen Species	Rickettsia conorii (strain ATCC VR-613 / Malish 7)
Purity	>85% (SDS-PAGE)
Sequence	MAELRLPPNS VVKKGREHKE QEEMLKPRKI KIIYRYDPDLN KNPTIDSFEI DLSKTGPMIL DALIKIKNEI DSTLTFRRSC REGICGSCAM NIDGTNTLAC IKPIEDISGD IKIYPLPHMK VVKDLVPDMS HFYAQYESIE PWLKTDSPTP SNSERLQSIK DREKLDGLYE CILCACCSTS CPSYWWNGDK YLGPAILLQA YRWIADSRDD NTGERLEALE DPFKLYRCHT IMNCTKTCPK GLNPAKAIGK IKSLIAERHG L
Source	Mammalian cell
Target Names	sdhB
Protein Names	Recommended name: Succinate dehydrogenase iron-sulfur subunit EC= 1.3.99.1
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.