



Recombinant Human N-acetyltransferase 6 (NAT6)

Product Code	CSB-BP821891HU
Abbreviation	NAT6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q93015
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MELILSTSPA ELTLDPACQP KLPLDSTCQP EMTFNPGPTE LTLDPHQPE ETPAPSLAEL TLEPVHRRPE LLDACADLIN DQWPRSRTSR LHSLGQSSDA FPLCLMLLSP HPTLEAAPVV VGHARLSRVL NQPQSLLVET VVVARALRGR GFGRRMEGL EVFARARGFR KLHLTTHDQV HFYTHLGYQL GEPVQGLVFT SRRLPATLLN AFPTAPSPRP PRKAPNLTAQ AAPRGPKGPP LPPPPPLPEC LTISPPVPSG PPSKSLLETQ YQNVGRPIF WMEKDI
Source	Baculovirus
Target Names	NAT6
Protein Names	Recommended name: N-acetyltransferase 6 EC= 2.3.1.- Alternative name(s): Protein fusion-2 Short name= Protein fus-2
Expression Region	1-286
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.