



Recombinant Syntrophin-1 (stn-1)

Product Code	CSB-BP821911CXY
Abbreviation	stn-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q93646
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MAAVRSGGLVD IFVQQQWHRV LATLDPTAIT LQTMEQNEAE AEKRTVRVVK YDGNGLGISI KGGRDNNMPI VISKIFKGMA ADQAGELFLD DVIISVNGEN LLDASHEEAV RALKRAGR VV DLQVQYRRED MMHRENIVEN VEWDDDIRER VRTIGLKLAY VARAGIDADA EGRILEMRSP SGRYSLAMRC SSSEEADGWF EALHACTTCL LTQALAQVNI MLGNNPQVRH MGWIAEQVSE NGISMWKPKF MTLTNSEILF YEAVPQLKAE WAEPRLVRPL VATRVVQTSS RTAPVIKGLT DVISFRMRTG TQQGVRHTI RVETHAELAR WVRVVIGGY EACLSTSQVS APCLWRGESC ELIVNLDNGI SLLSSTGEVL WQHSFETIRA TGDDGGRFLW VDFGPPHGEQ ELDLLNSAKP VVFILHSFLA TKVYRLGLYA
Source	Baculovirus
Target Names	stn-1
Protein Names	Recommended name: Syntrophin-1
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.